


A Personal Sampler to Estimate Particle Deposition in the Human Respiratory Tract

Kirsten Koehler, Phillip Clark, John Volckens

Department of Environmental and Radiological Health Sciences, Colorado State University, Fort Collins, CO

Introduction



*Measurement of aerosol deposition is a better surrogate for dose, risk assessment.

Advantages of Estimating Deposition vs. Exposure

- *Not all inhaled particles deposit in the respiratory tract.
- *Filter based samplers measure all particles (inhaled and exhaled).
- *This results in error
 - *Deposition vs. Inhalability
 - *Deposition vs. PM_{2.5}
- *A sampler that mimics lung deposition may minimize this error.




Figure 1: Schematic showing error between an inhalable sampler and a deposition sampler for a typical aerosol size distribution. The inhalable sampler captures all the particles, including those exhaled.

Figures 1 and 2 exemplify the difficulty in relating inhalability (or PM exposure) with deposition, dose, and risk.

Measurement Error: Exposure vs. Deposition





Figure 2: Error between common filter-based samplers and total respiratory deposition. (a) IOM sampler. Ovals represent particle size ranges reported in the literature for occupational exposures. (b) PM_{2.5} sampler.

Aerosol Deposition Model

Lung deposition sampler engineered from porous foam.

$$D = 1 - \exp \left\{ - \frac{t}{d_f} \left[a(St)^b + c(Ng)^d + e(Pe)^{-f} \right] \right\} \quad \text{Clark et al. (2009)}$$

Foam Properties Impaction Settling Diffusion





Figure 3: Particle deposition fraction for the foam sampler compared to ICRP (1994) total deposition and inhalable fraction. The red line represents modeled data; the symbols represent experimental data. Error bars represent one standard deviation.

*We inverted the model to find foam parameters that give the best agreement to deposition fractions given by the International Commission on Radiological Protection (ICRP model, 1994).

*We have verified that loadings of up to 20 mg have no effect on deposition efficiency.

References:

ICRP. (1994) Human Respiratory Tract Model for Radiological Protection. Oxford, UK:Elsevier Science, Ltd.
Clark, P. et al. An improved model for particle deposition in porous foams, J. Aerosol Sci., 2009; 40(7): 563-572.

Designing The Foam Sampler

Design Guidelines

	Pump	Foam	Battery	Total
Inexpensive	\$50	\$1	\$30	\$81
Lightweight	150 g	0.05 g	185 g	335 g
Reliable	Figure 6	Pressure drop constant with aerosol loading	Figure 6	

Using off-the-shelf parts the deposition sampler was built for < \$100 and weighs about 0.3 kg (0.5 lb).




Figure 4: Pump, battery and control electronics of the sampling system. The entire unit is small enough to fit in a fanny pack.



Figure 5: The sampler can be worn on the shoulder near a workers breathing zone.

Battery Operated Pump Performance

- *Tested variation in flow over a typical 8-hr exposure period.
- *Linear decrease in flow rate with running time, due to the battery.
- *The ~10% decrease in flow has a negligible effect on the deposition efficiency of the foam.




Figure 6: Pump flow rate and battery voltage (after potentiometer) as a function of running time of the battery.

Known Limitations

- *Gravimetric analysis influenced by humidity, electrostatics.

Conclusions

- *A foam sampling system was developed that mimics the aerosol deposition in the human respiratory tract and is expected to provide a more physiologically relevant estimation of risk.
- *This sampler presents a significant reduction in cost because inexpensive pumps can be used to control flow.
- *Some difference in total and deposited aerosol concentration are noted for common exposures.

Acknowledgements

This work was funded by grant #OH009248 from the National Institute for Occupational Safety and Health and by Grant Number OH009229 from CDC NIOSH Mountain and Plains Education and Research Center

Real-Time Exposure Monitoring

- *Can apply on a real-time basis using a condensation particle counter or nephelometer (here, Personal DataRam, pDR, Thermo Scientific, Inc.).
- *Unmodified pDR gives the total aerosol concentration in mg/m³.
- *pDR with the foam insert measures the exhaled portion.
- *Deposited Aerosol = Total Aerosol – Exhaled Aerosol.




Figure 7: Time series of Total and Deposited aerosol concentration for several different aerosol exposures.

Contact Information:

Kirsten.Koehler@colostate.edu John.Volckens@colostate.edu

More details on this work can be found in Koehler et al. (2009):

Development of a Sampler for Total Aerosol Deposition in the Human Respiratory Tract, Ann. Hyg., 2009; 53: 731-738.